Multi-channel Pumps

More efficiency, even more possibilities

With the easily exchangeable cassettes, the throughput of the Hei-FLOW multi-channel pump can be increased to up to 12 simultaneously operated channels.

The following models are suitable for multi-channel operation: **Hei-FLOW Core/Expert/Ultimate 600**

Hei-FLOW UN

Simply select the adapter and multi-channel pump head for the appropriate Hei-FLOW model and fit them with suitable cassettes and tubings.


- When using tubes with different diameters per cassette, up to 12 individual pump volumes can be processed in one working part
- The tube can be changed easily and in a matter of seconds
- Pump heads with 8-roller system are also available to reduce pulsation
- A snap-action device makes inserting all cassettes child's play and even allows easy replacement during ongoing operation

Flow rates of individual tubing sizes for multi-channel pump heads

Multi-channel pump head C 4

For Cassette small








Multi-channel pump head C 12

For Cassette small






Multi-channel pump head C 8

For medium or large cassettes





Multi-channel pump heads

Easy to configure or retrofit

Multi-channel cassettes

The brackets of the 2-stopper tubings prevent the slipping of the tube when using the small cassette.



Easily exchangeable cassettes even during the pumping process. The roller contact pressure is adjusted by means of an adjusting screw. Different tubing and sizes can be used in each cassette.

Precise dosing or customized pumping

Low-pulsation pumping with the C 4 and C 12 multi-channel pump heads thanks to the 8-roller system and high-precision dosing depending on the tubing configuration. The C 12 model is optimally equipped for the smallest volumes thanks to an integrated gear support – for flow rates from 0.005–54 ml/min. For Cassette small (C 4/C 12), Two-Stop tubing is required. For Cassette medium and Cassette large (C 8) tubing by the meter.



Multi-channelpump head C 4

- Can be equipped with 4 x Cassette small
- 8 rollers for low-pulsation pumping

P/N 524-80420-00



Multi-channelpump head C 8

- Can be equipped with 8x Cassette medium or 4 x Cassette large
- 4-roller system

P/N 524-40810-00



Multi-channelpump head C 12

- Can be equipped with 12 x Cassette small
- Due to integrated gear reduction ideal for pumping smallest volumes
- 8 rollers for low-pulsation pumping

P/N 524-81220-00



Flow rates from

0.8 and 1.7 mm

meter

Cassette small

- Flow rates from 0.005-37.0 ml/min
- Suitable for tubes with 0.9 mm tube wall thickness
- Available tube diameters: 0.2/0.5/0.9/1.4 and 2.8 mm
- Special tube piece with 2 stoppers (length 40 cm) required for insertion into the cassette
- The tube is fixed by tube stoppers
- With tubing connectors and extension tubings, it is possible to extend the tubing length by the meter

Equipped with:

Equipped with:

Multi-channel pump head C 4: max. 4 x Cassette small

Multi-channel pump head C 12: max. 12 x Cassette small

P/N 524-90022-00

Multi-channel pump head C 8: max. 8x Cassette medium

P/N 524-90021-00



Cassette medium

0.22-25.0 ml/min Suitable for tubes with 1.6 mm tube wall thickness

• Available tube diameters:

Tubes available by the



Cassette large

- Flow rates from 1.0-277.0 ml/min
- Suitable for tubes with 1.6 mm tube wall thickness
- Available tube diameters: 1.7/3.1/4.8 and 6.4 mm
- Tubes available by the meter

Equipped with:

Multi-channel pump head C 8: max. 4 x Cassette large

P/N 524-90010-00

Tubing sizes for multi-channel pump heads

Tube sizes		•	•	0	0	0
Inside diameter	mm	0.25	0.51	0.89	1.42	2.79
Outside diameter	mm	2.05	2.31	2.69	3.22	4.59
Tube wall thickness (wt)	mm	0.9	0.9	0.9	0.9	0.9
Max. operating pressure (duration/short-term)	bar	0.5/1.5	0.5/1.5	0.5/1.5	0.5/1.5	0.5/1.5
Suction lift	mH ₂ O	7	7	7	7	7

Tube sizes		0	0	0	0	0
Inside diameter	mm	0.8	1.7	3.1	4.8	6.4
Outside diameter	mm	4	4.9	6.3	8	9.5
Tube wall thickness (wt)	mm	1.6	1.6	1.6	1.6	1.6
Max. operating pressure (duration/short-term)	bar	0.7/1.7	0.7/1.7	0.7/1.7	0.7/1.7	0.5/1.5
Suction lift	mH ₂ O	8.8	8.8	8.8	8.8	6.7

Mean value of the flow rate in combination with pump head and pump drive

Hei-FLOW Core 120 Hei-FLOW Expert 120 Hei-FLOW Ultimate 120		min.	max.	min.	max.	min.	max.	min.	max.	min.	max.	max. number of cass.
Cassette small Pump head C 12	ml/min	0.005	0.11	0.01	0.54	0.03	1	0.10	3	0.29	9	12
Cassette small Pump head C 4	ml/min	0.02	0.49	0.08	2	0.24	6	0.60	14	2	36	4

All flow rate data refer to $\mathsf{Tygon}^{\texttt{B}}$ (standard) tubes and the medium water.

Order numbers

Tubing	P/N	P/N	P/N	P/N	P/N
Silicone					
Two-Stop tubing for Cassette small			525-30014-00	525-30015-00	525-30016-00
Extension tubes (by the meter)			525-30024-00	525-30025-00	525-30026-00
Viton®					
Two-Stop tubing for Cassette small			525-00014-00	525-00015-00	525-50016-00
Extension tubes (by the meter)			525-00024-00	525-00025-00	525-50026-00
PharMed®					
Two-Stop tubing for Cassette small	525-20012-00	525-20013-00	525-20014-00	525-20015-00	525-20016-00
Extension tube (by the meter)	525-20022-00	525-20023-00	525-20024-00	525-20025-00	525-20026-00
Tygon [®] standard					
Two-Stop tubing for Cassette small	525-60012-00	525-60013-00	525-60014-00	525-60015-00	525-60016-00
Extension tubes (by the meter)	525-60022-00	525-60023-00	525-60024-00	525-60025-00	525-60026-00
Tube connector (PTFE)	526-22000-00	526-22000-00	526-22000-00	526-22000-00	526-22000-00

Hei-FLOW Core 120 Hei-FLOW Expert 120 Hei-FLOW Ultimate 120		min.	max.	max. number of cass.								
Cassette medium Pump head C 8	ml/min	0.22	6.8	1	25							8
Cassette large Pump head C 8	ml/min			1	25	3.7	88.5	7.7	184	11	277	4

All flow rate data refer to $\mathsf{Tygon}^{\textcircled{B}}$ (standard) tubes and the medium water.

Tubing	P/N	P/N	P/N	P/N	P/N
Silicone	525-33000-00	525-34000-00	525-36000-00	525-30027-00	525-30028-00
Viton®	525-53000-00	525-54000-00	525-56000-00	525-50027-00	525-50028-00
PharMed [®]	525-23000-00	525-24000-00	525-26000-00	525-20027-00	525-20028-00
Tygon [®] standard	525-63000-00	525-64000-00	525-66000-00	525-60027-00	525-60028-00
Tygon [®] for hydrocarbons	525-73000-00	525-74000-00	525-76000-00	525-70027-00	525-70028-00
Tygon [®] 2001 for food	525-83000-00	525-84000-00	525-86000-00	525-80027-00	525-80028-00