


InLab® 741

Cond °C ISM



For low conductivity range samples
like pure water and ultra pure water.



Sensor type	conductivity cell
Order number	51344024 - cable length 1.2 m 51344026 - cable length 5 m
Specifications	
measuring range	0.001...500 $\mu\text{S}/\text{cm}$
zero point (pH)	-
temperature range	0...100 °C
temperature probe	NTC 30 k Ω
shaft material	stainless steel
membrane glass	-
membrane resistance	-
reference system	-
reference diaphragm	-
reference electrolyte	-
bridge electrolyte	-
electrode head	-
connector	Mini-DIN
cell constant (cm^{-1})	0.105 (certified)
cell type	2 steel poles
indicating metal	-
Storage	
dry	