


# InLab® 742-ISM

Cond °C ISM



For portable applications in the lab and outdoor. Ideal for low conductivity range samples like pure water and ultra pure water.



<b>Sensor type</b>	conductivity cell
<b>Order number</b>	51344116 - cable length 1.8 m 51344118 - cable length 5 m
<b>Specifications</b>	
measuring range	0.001...500 µS/cm
zero point (pH)	-
temperature range	0...100 °C
temperature probe	NTC 30 kΩ
shaft material	stainless steel
membrane glass	-
membrane resistance	-
reference system	-
reference diaphragm	-
reference electrolyte	-
bridge electrolyte	-
electrode head	-
connector	LTW (IP67)
cell constant (cm <sup>-1</sup> )	0.105 (certified)
cell type	2 steel poles
indicating metal	-
<b>Storage</b>	
dry	