


InLab[®] Redox Flow

Combined ORP electrode with platinum ring

Scale 1:1



Sensor type combined ORP electrode

Order number 51343201

Specifications

pH range -

zero point (pH) -

temperature range 0...100 °C

temperature probe -

shaft material glass

membrane glass -

membrane resistance -

reference system ARGENTHAL[™]

reference diaphragm movable glass sleeve

reference electrolyte 3 mol/L KCl

bridge electrolyte -

electrode head S7

cable length -

connector -

cell constant (cm⁻¹) -

cell type -

indicating metal platinum ring

Storage

InLab storage solution (30111142)

or 3 mol/L KCl (51350072)