


# InLab<sup>®</sup> Redox-L

Long ORP electrode with platinum ring

Scale 1:1



**Sensor type** combined ORP electrode

**Order number** 51343202

**Specifications**

pH range -  
zero point (pH) -

temperature range 0...100 °C  
temperature probe -

shaft material glass  
membrane glass -  
membrane resistance -

reference system ARGENTHAL™  
reference diaphragm ceramic  
reference electrolyte 3 mol/L KCl  
bridge electrolyte -

electrode head S7  
cable length -  
connector -

cell constant (cm<sup>-1</sup>) -  
cell type -

indicating metal platinum ring

**Storage**

InLab storage solution (30111142)  
or 3 mol/L KCl (51350072)