

Accessories



3005



3006



Shaking tray

made of stainless steel, 330 x 330 mm, with holes to accommodate clamps for Erlenmeyer flasks and other accessory equipment

Order No. 3950




Clamps

made of stainless steel for Erlenmeyer flasks, for screwing on to shaking tray 3950, complete with attachment materials

- Order No. 3983** for 25 ml flasks (45*)
- Order No. 3984** for 50 ml flasks (25*)
- Order No. 3985** for 100 ml flasks (16*)
- Order No. 3986** for 200 ml flasks (12*)
- Order No. 3987** for 250-300 ml flasks (9*)
- Order No. 3988** for 500 ml flasks (9*)
- Order No. 3989** for 1000 ml flasks (4*)
- Order No. 3990** for 2000 ml flasks (2*)

Other clamps on request

* max. number of clamps per shaking tray




Test tube rack

made of stainless steel, for test tubes; perforated section can be tilted by 90°, with clamping springs for secure grip and silent shaking of glasses; can be screwed on to shaking tray 3950

Order No. 3953

for a max. of 24 tubes: 12 - 17 mm diameter
Capacity of tray: three test tube racks



Clamps

made of stainless steel, for separating funnels, for screwing on to shaking tray 3950

Order No. 3957 for 50 ml separating funnels (6*)
(ISO/Squibb, neck 19/26)

Order No. 3958 for 100 ml separating funnels (6*)
(ISO/Squibb, neck 19/26)


Order No. 3959 for 250 ml separating funnels (4*)
(ISO/Squibb, neck 19/26)

Order No. 3955 for 250 ml separating funnels (4*)
(conical form, neck 29/32)

Order No. 3956 for 500 ml separating funnels (3*)
(ISO/Squibb, neck 29/32)

Other clamps on request

* Max. number of clamps per shaking tray




Holding frame

made of stainless steel, for attachment of a test plate, for screwing on to shaking tray 3950, complete with attachment materials

Order No. 3910

Capacity of tray: four test plates



Non-slip mat

for the shaking platform, 300 x 300 mm size, for slow moving, e.g. nutrient solutions in Petri dishes and flasks

Order No. 3951



Universal mount

for secure attachment of different shaking objects between the four rubber-coated bars

Order No. 3952

Accessories



3011




3012



3013




3014



Order No. 3965

Non-slip mat


for the shaking platform, 420 x 420 mm size, for slow moving, e.g. nutrient solutions in Petri dishes



Order No. 3967

Universal mount

for secure attachment of different shaking objects between the six rubber-coated bars



Order No. 3966

Shaking tray

made of stainless steel, 450 x 450 mm, with holes to accommodate clamps for Erlenmeyer flasks and other accessory equipment



Order No. 3968

Platform frame

with four levels, with top rims, three levels removable. Frame and trays made of stainless steel. Distance between trays - 60 mm. Grooved rubber mats on trays to ensure secure grip of objects.



Clamps

made of stainless steel for Erlenmeyer flasks, for screwing on to shaking tray 3966, complete with attachment materials

- Order No. 3983** for 25 ml flasks (79*)
- Order No. 3984** for 50 ml flasks (49*)
- Order No. 3985** for 100 ml flasks (36*)
- Order No. 3986** for 200 ml flasks (22*)
- Order No. 3987** for 250-300 ml flasks (16*)
- Order No. 3988** for 500 ml flasks (12*)
- Order No. 3989** for 1000 ml flasks (9*)
- Order No. 3990** for 2000 ml flasks (4*)

Other clamps on request

* Max. number of clamps per shaking tray

Accessories



3015




3017



3016



3018



Shaking tray

made of stainless steel, 450 x 450 mm, with holes to accommodate clamps for Erlenmeyer flasks and other accessory equipment

Order No. 3966




Clamps

made of stainless steel for Erlenmeyer flasks, for screwing on to shaking tray 3966, complete with attachment materials

- Order No. 3983** for 25 ml flasks (79*)
- Order No. 3984** for 50 ml flasks (49*)
- Order No. 3985** for 100 ml flasks (36*)
- Order No. 3986** for 200 ml flasks (22*)
- Order No. 3987** for 250-300 ml flasks (16*)
- Order No. 3988** for 500 ml flasks (12*)
- Order No. 3989** for 1000 ml flasks (9*)
- Order No. 3990** for 2000 ml flasks (4*)

Other clamps on request

* Max. number of clamps per shaking tray



Test tube rack

made of stainless steel, for test tubes; perforated section can be tilted by 90°, with clamping springs for secure grip and silent shaking of glasses; can be screwed on to shaking tray 3966

Order No. 3953

for a max. of 24 tubes: 12 - 17 mm diameter
capacity of tray: six test tube racks



Clamps

made of stainless steel, for separating funnels, for screwing on to shaking tray 3966

Order No. 3957 for 50 ml separating funnels (11*)
(ISO/Squibb, neck 19/26)

Order No. 3958 for 100 ml separating funnels (11*)
(ISO/Squibb, neck 19/26)

Order No. 3959 for 250 ml separating funnels (8*)
(ISO/Squibb, neck 19/26)

Order No. 3955 for 250 ml separating funnels (8*)
(conical form, neck 29/32)

Order No. 3956 for 500 ml separating funnels (6*)
(ISO/Squibb, neck 29/32)

Other clamps on request

* Max. number of clamps per shaking tray


Holding Frame

made of stainless steel, for attachment of a test plate, for screwing on to shaking tray 3966, complete with attachment materials

pict. see Shakers 3019 - 3020

Order No. 3910

Capacity of tray: six test plates



Non-slip mat

for the shaking platform, 420 x 420 mm size, for slow moving, e.g. nutrient solutions in Petri dishes

Order No. 3965




Universal mount


for secure attachment of different shaking objects between the six rubber-coated bars

Order No. 3967

Accessories



3019



3020



Order No. 3980

Shaking tray

made of 3-mm-thick anodised aluminium, 670 x 537 mm, with holes to accommodate clamps for Erlenmeyer flasks and other accessory equipment




Clamps

made of stainless steel for Erlenmeyer flasks, for screwing on to shaking tray 3980, complete with attachment materials

- Order No. 3983** for 25 ml flasks (99*)
- Order No. 3984** for 50 ml flasks (99*)
- Order No. 3985** for 100 ml flasks (50*)
- Order No. 3986** for 200 ml flasks (26*)
- Order No. 3987** for 250-300 ml flasks (26*)
- Order No. 3988** for 500 ml flasks (26*)
- Order No. 3989** for 1000 ml flasks (12*)
- Order No. 3990** for 2000 ml flasks (9*)

Other clamps on request

* Max. number of clamps per shaking tray



Order No. 3953

for a max. of 24 tubes: 12 - 17 mm diameter
Capacity of tray: nine test tube racks

Test tube rack

made of stainless steel for test tubes; perforated section can be tilted by 90°, with clamping springs for secure grip and silent shaking of glasses; can be screwed on to shaking tray 3980



Order No. 3910

Capacity of tray: fifteen test plates

Holding frame

made of stainless steel, for attachment of a test plate, for screwing on to shaking tray 3980, complete with attachment materials



Order No. 3982

Universal mount

for secure attachment of shaking objects. The mount consists of two longitudinal bars and six transverse bars. Screw-type clamps are used for adjusting the longitudinal bars vertically and transverse bars horizontally.



Order No. 3981

Rack frame

for doubling the capacity by accommodating two shaking trays 3980; distance between trays - 205 mm. The bottom tray can be loaded with flasks up to a max. of 500 ml.