


Heating mantle | beaker

for beaker, integrated controller, up to 450°C



WHM for beaker, 450°C integrated controller

Ideal for:


- for beaker 100 - 5.000 ml

Features:

- 100°C - 450°C
- full choice of vessel, capacity, housing and control are offered
- nickel chrome heating element
- coil type: acid-resistant, long life time
- all kinds of beaker can be used (Duran®/Kimax/Pyrex/Iwaki)
- CE certified and unique serial number for tracing

Controller:

- built-in *temperature controller*



1.000 ml beaker heating mantle

| Application (beaker brand) | Beaker capacity (ml) | Beaker Ø (mm) | Watts (W) | Inside depth (mm) | Mantle outside Ø x H mm | Support rod Ø x H mm | Order number |
|-------------------------------|----------------------------|------------------|--------------|-------------------------|-------------------------------|-------------------------|--------------|
| | | | | | | | with 230V |
| common-use | 100 | 50 - 52 | 80 | 40 | Ø 145 x 105 | Ø 13 x 440 | DH.WHM12141 |
| common-use | 250 | 68 - 70 | 140 | 65 | Ø 145 x 140 | Ø 13 x 440 | DH.WHM12142 |
| for DU. | 400 | 79 - 81 | 240 | 70 | Ø 145 x 140 | Ø 13 x 440 | DH.WHM12143 |
| for Iw./Ki./Py. | 400 | 76 - 78 | 240 | 75 | Ø 145 x 140 | Ø 13 x 440 | DH.WHM121431 |
| common-use | 600 | 88 - 90 | 240 | 75 | Ø 145 x 140 | Ø 13 x 440 | DH.WHM12144 |
| common-use | 800 | 98 - 100 | 300 | 90 | Ø 190 x 165 | Ø 13 x 440 | DH.WHM12145 |
| common-use | 1,000 | 105 - 108 | 350 | 90 | Ø 190 x 165 | Ø 13 x 440 | DH.WHM12146 |
| common-use | 2,000 | 130 - 132 | 450 | 125 | Ø 190 x 180 | Ø 13 x 660 | DH.WHM12147 |
| for DU./Iw. | 3,000 | 151 - 153 | 550 | 160 | Ø 270 x 235 | Ø 13 x 660 | DH.WHM12148 |
| for Ki./Py. | 4,000 | 160 - 162 | 700 | 160 | Ø 270 x 235 | Ø 13 x 660 | DH.WHM121481 |
| for DU. | 5,000 | 169 - 171 | 800 | 160 | Ø 270 x 235 | Ø 13 x 660 | DH.WHM12149 |

Remarks: DU. = Duran®, Ki. = Kimax, Iw. = Iwaki, Py. = Pyrex