

Vacuum pump (oil-free)

ROCKER
make filtration easy

Rocker 800 / 810 up to -735mmHg, 100 l/min, 68dB



Rocker 800



Rocker 810

Ideal for:

- vacuum drying
- vacuum oven
- high volume degas
- high volume vacuum filtration

Features:

- lubricant and oil free piston
- no air pollution
- maintenance-free
- Prote™ filter prevents from liquid being sucked into the pump and filters particles as well as moisture
- CE certified

Safety mechanism:

- thermal protection device ensures shutting down the pump when overheated and resumes working when the temperature cooled down

Model	Rocker 800		Rocker 810	
Max. power (W)	220	200	220	200
Max. current (A)	2.0	0.9	2.0	0.9
Max. vacuum	-670 mmHg (119 mbar)		-735 mmHg (33 mbar)	
Max. flow rate	100 l/min	84 l/min	60 l/min	50 l/min
Moisture trap / vacuum regulator	Yes		Yes	
Overflow protection	No		No	
Motor ration	1750 rpm	1450 rpm	1750 rpm	1450 rpm
Horse power	1/3 HP		1/3 HP	
Noise level	68 dB		68 dB	
Port thread	3/8 inch (10 mm)		3/8 inch (10 mm)	
Net weight	7.8 kg		7.8 kg	
Dimension (LxWxH)	34.5 x 17.5 x 20 cm		34.5 x 17.5 x 20 cm	
Power supply	1 Phase AC 120V, 60 Hz or AC 230V, 50/60 Hz			
Order number 110V	00829 167800-11	-	00829 167810-11	-
Order number 230V	-	00829 167800-22	-	00829 167810-22