

Vacuum pump (oil-free)

ROCKER
make filtration easy

Rocker 801 / 811

up to -735 mmHg, 80 l/min, 68dB



Rocker 801



Rocker 811

Ideal for:

- filtration for microbiology test
- air sampling
- suspended solid test
- vacuum filtration

Features:

- lubricant and oil free piston
- no air pollution
- maintenance-free
- Prote™ filter prevents from liquid being sucked into the pump and filters particles as well as moisture
- CE certified

Safety mechanism:

- thermal protection device ensures shutting down the pump when overheated and resumes working when the temperature cooled down

| Model | Rocker 801 | | Rocker 811 | |
|----------------------------------|---|-----------------|---------------------|-----------------|
| Max. power (W) | 220 | 200 | 220 | 200 |
| Max. current (A) | 2.0 | 0.9 | 2.0 | 0.9 |
| Max. vacuum | -670 mmHg (119 mbar) | | -735 mmHg (33 mbar) | |
| Max. flow rate | 80 l/min | 68 l/min | 48 l/min | 40 l/min |
| Moisture trap / vacuum regulator | Yes | | Yes | |
| Overflow protection | No | | No | |
| Motor ration | 1750 rpm | 1450 rpm | 1750 rpm | 1450 rpm |
| Horse power | 1/3 HP | | 1/3 HP | |
| Noise level | 68 dB | | 68 dB | |
| Port thread | 3/8 inch (10 mm) | | 3/8 inch (10 mm) | |
| Net weight | 7.9 kg | | 7.9 kg | |
| Dimension (LxWxH) | 34.5 x 17.5 x 20 cm | | 34.5 x 17.5 x 20 cm | |
| Power supply | 1 Phase AC 120V, 60 Hz or AC 230V, 50/60 Hz | | | |
| Order number 110V | 00829 167801-11 | - | 00829 167811-11 | - |
| Order number 230V | - | 00829 167801-22 | - | 00829 167811-22 |