


# InLab<sup>®</sup> Reach 225

Long sensor for general applications

pH Ref °C ISM



Scale 1:1



**Sensor type** combined pH electrode

**Order number** 30244733

## Specifications

pH range 0...14

zero point (pH) 7.00 ± 0.25

temperature range 0...100 °C

temperature probe -

shaft material glass

membrane glass HA

membrane resistance < 600 MΩ

reference system ARGENTHAL™

reference diaphragm ceramic

reference electrolyte 3 mol/L KCl

bridge electrolyte -

electrode head S7

cable length -

connector -

cell constant (cm<sup>-1</sup>) -

cell type -

indicating metal -

## Storage

InLab storage solution (30111142)

or 3 mol/L KCl (51350072)