

Filter holder | PES

Filter Ø47mm, 300ml / 500ml magnetic PES filter holder

Features:

- magnetic connection facilitates fast and stable installation without clamp
- strong magnetic seal allows one-handed operation, prevents possible twisting and tearing of the membrane
- membrane guide on the base helps positioning membrane easily
- PES structure offers excellent smash and high temperature resistance (up to 180°C)
- autoclavable
- reusable
- funnel graduated in 50ml increments for easy measurement



MF 3 - 300ml
Ø47mm, with stopper,
long stem




MF 5 - 500ml
Ø47mm, with stopper, long stem



MF 3a - 300ml
Ø47mm, short stem



MF 5a - 500ml
Ø47mm, short stem




MF 3 / 5

MF 3 / 5 magnetic PES filter holder, stopper, long stem

#	Item	Material	Size	Order number
	300ml, 47mm magnetic PES filter holder, complete			00829 200300-01
	500ml, 47mm magnetic PES filter holder, complete			00829 200500-01
1	Lid kit	PP		00829 200300-36
2	Funnel	PES	300 ml 500 ml	00829 200300-30 00829 200500-30
3	Support base	PES	long stem	00829 200300-32
4	Stopper	Silicone	No. 8	00829 167110-16

MF 3a / 5a magnetic PES filter holder, short stem

#	Item	Material	Size	Order number
	300ml, 47mm magnetic PES filter holder, complete			00829 200310-01
	500ml, 47mm magnetic PES filter holder, complete			00829 200510-01
1	Lid kit	PP		00829 200300-36
2	Funnel	PES	300 ml 500 ml	00829 200300-30 00829 200500-30
3	Support base	PP	short stem	00829 200310-32



MF 3a / 5a