

Pump accessories

1200 ml PES waste bottle

Specification	
Port thread (I.D. tubing)	5/16 inch (8 mm)
Graduated increment	100 ml
Overflow protection	Yes
Cap with suction fitting	Yes
Material	Container: PES Cover: PES Overflow prot.: PP
Order number	00829 197000-11



1200 ml PC waste bottle

Specification	
Port thread (I.D. tubing)	5/16 inch (8 mm)
Graduated increment	100 ml
Overflow protection	Yes
Cap with suction fitting	Yes
Material	Container: PC Cover: PES Overflow prot.: PP
Order number	00829 197000-11-PC



3000 ml PC waste bottle with grip

Specification	
Port thread (I.D. tubing)	5/16 inch (8 mm)
Graduated increment	200 ml
Overflow protection	Yes
Material	Container: PC Overflow prot.: PP
Order number	00829 167240-33



1200 ml PES waste bottle

Specification	
Port thread (I.D. tubing)	5/16 inch (8 mm)
Graduated increment	100 ml
Overflow protection	Yes
Material	Container & cover: PES Overflow prot.: PP
Order number	00829 197000-12



1000 ml / 3000 ml PC waste bottle

Specification	1000 ml	3000 ml
Port thread (I.D. tubing)	5/16 inch (8 mm)	
Graduated increment	100 ml	
Overflow protection	Yes	
Material	Container: PC Overflow prot.: PP	
Order number	00829 167200-31	00829 167200-33



2000 ml / 4000 ml HPDE waste bottle

Specification	2000 ml	4000 ml
Port thread (I.D. tubing)	1/4 inch (6.5 mm)	
Graduated increment	500 ml	
Overflow protection	No	Yes
Material	Cover: PP	
Order number	00829 167200-51	00829 167200-46



2000 ml / 4000 ml PP waste bottle

Specification	2000 ml	4000 ml
Port thread (I.D. tubing)	1/4 inch (6.5 mm)	
Graduated increment	500 ml	
Overflow protection	No	Yes
Material	Container: PP Cover: PP Overflow prot.: PP	
Order number	00829 167200-50	00829 167200-47



4000 ml PP waste bottle with drain

Specification	
Port thread (I.D. tubing)	1/4 inch (6.5 mm)
Graduated increment	500 ml
Overflow protection	Yes
Drain	Yes
Material	Container: PP Cover: PP Overflow prot.: PP Tube: rubber
Order number	00829 167200-48



1000 ml receiver flask with drain & silicone fixing sucker

Specification	
Port thread (I.D. tubing)	5/16 inch (8 mm)
Graduated increment	100 ml
Drain	Yes
Silicone fixing sucker	Yes
Stopper	No. 8 (optional)
Material	Container: Borosilicate glass Connector: PP
Order number	00829 167100-15



Fixing sucker for receiving bottles

Specification	
Suitable flask size	1000 ml
Material	Silicone
Order number	00829 167200-45

