

Magnetic stirrer with hotplate | Hotplate

Accessories: external temperature probe, stand rods, holder/clamp/clip

Temperature probes for hotplates and magnetic stirrer with hotplate

For device(s)	Design	Measurements	Article	Order number
HP-500D MSH-500D MSH-50D HP-20D HP-30D MSH-20D MSH-30D SMHS	Standard	Ø4xL250 mm, cable: 190 mm	SS500	DH.WHP503060
HP-500D MSH-500D MSH-50D HP-20D HP-30D MSH-20D MSH-30D SMHS	PTFE-coated	Ø4xL250 mm, cable: 190 mm	SS510	DH.WHP503061



Stand rods for hotplates and magnetic stirrer with hotplate

For device(s)	Design	Measurements	Article	Order number
HP-30A, HP-30D MSH-30A, MSH-30D	without thread	Ø12.7xh450 mm	RD100	DH.WHP503011
HP-20A, HP-20D, HP-500D MSH-20A, MSH-20D MSH-50A, MSH-50D MSH-500D, SMHS	with thread	Ø12xh450 mm	RD200	DH.WHP503036



Holder, clamp and clip suitable for stand rods RD100 and RD200

For device(s)	Article	Order number
all hotplates and magnetic stirrer with hotplate	CL220	DH.WHP503012

