QualiRod cylindrical





Cylindrical hollow drill with a sharp cutting edge, can easily be driven into paste-like and semi-solid substances. The samples cut out are subsequently removed with the appropriate push-rod ejector. Total length 260 mm. Inner \emptyset 8 mm, outer \emptyset 10 mm.

- ▶ Stainless steel AISI 304 (1.4301)
- Sterilisable
- ▶ Inner Ø 8 mm, outer Ø 10 mm

QualiRod cylindrical Drilling depth mm	Quantity scale from	Item no.
200 Accessories	5	5311-0101 N4
Description		Item no.
Transport case, inner dimensions 258x178x45 mm		3600-0021 N3