









SteriBag StandUp is a sterile sample bag with a stand-up base. The special bottom unfolds as soon as the bag is filled. This makes SteriBag StandUp stable and able to stand on its own without a separate stand or support, leaving both your hands free to take the sample and fill the bag. The wire bracket at the closure keeps the bag open, making it easy and convenient to fill.

The sample bag is the ideal space-saving, unbreakable, inexpensive and sterile replacement for bottles and beakers.

Can be used for samples available as liquids, powders, solids and pastes; suitable for foods and can be frozen.

The stable standing SteriBag StandUp sample bag is easy and safe to transport.

SteriBag StandUp, the sample bag for easy handling of samples.

- Stands stably
- Protected against discharge
- ▶ Sterile, RNase, DNase and pyrogen-free
- With labelling field

SteriBag StandUp sample bag			PA PE	sterilised sterilised
Capacity ml	L x W mm	Pieces per pack		Item no.
400	229 x 114	500		5344-0400 N3
600	229 x 152	500		5344-0600 N3