


Crucible tongs with ridges



Crucible tong with ridges, stainless steel



Material	A (Length)	a (Length)	d/Ø	D/Ø	Article no.
18/10 E-Poli	160	85	18	4	2099
18/10 E-Poli	200	95	20	4,5	2100
18/10 E-Poli	220	100	20	4,5	2101
18/10 E-Poli	250	115	30	4,5	2102
18/10 E-Poli	300	145	30	5	2103
18/10 E-Poli	400	175	35	6	2104
18/10 E-Poli	500	225	35	6	2105
18/10 E-Poli	600	270	40	8	2106
18/10 Stainless	200	95	20	4,5	2130
18/10 Stainless	220	100	20	4,5	2131
18/10 Stainless	250	115	30	4,5	2132
18/10 Stainless	300	145	30	5	2133
PTFE coating	200	95	20	4,5	2400
PTFE coating	250	115	30	4,5	2401