

PRODUCT DATA SHEET

VS 15 R

The robust low-noise rocking shakers from the Varioshake product line are sure to impress with a reduced footprint in the laboratory as well as load-independent speed and digital control. They are perfect for continuous operation.

Technical Attributes	VS 15 R
Ambient temperature min.	10 °C
Ambient temperature max.	50 °C
Shaking platform width	450 mm
Shaking platform depth	450 mm
Load capacity max.	15 kg
Shaking amplitude	3 degrees
Shaking frequency min.	2 min-1
Shaking frequency max.	50 min-1
Shaking motion	rocking
Timer	1 minute - 99:59 hours or continuous operation
Dimensions (WxDxH) in mm	480x487x160
Weight	19.5 kg

LAUDA DR. R. WOBSER GMBH & CO. KG Pfarrstraße 41/43 • 97922 Lauda-Königshofen Postfach 1251 • 97912 Lauda-Königshofen • DE

T +49 (0) 9343 503-0 · F +49 (0) 9343 503-222 info@lauda.de · www.lauda.de WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen Registergericht Mannheim • HRA 560069

Persönlich hattende Gesellschatterin: LAUDA DR. R. WOBSER Verwaltungs-GmbH Sitz Lauda-Königshofen Registergericht Mannheim • HRB 560226 Geschäftsführer: Dr. Gunther Wobser (Vors.), Dr. Mario Englert Dr. Marc Stricker Beirat: Dr. Gerhard Wobser