The professional way for your pipetting applications

WITOPET professional single-channel is the classical microliter pipette, made for the daily use in the laboratory: accurate, robust and handy. The seven models available cover a volume range from 0.2 µl to 5 ml.



Pipette support

suitable for up to seven WITOPET professional pipettes 5 406 107

WITOPET professional



Variable single-channel

Your advantage at a glance:

- · high quality pipette, made in Germany
- · economical priced
- perfect for precise laboratory applications
- 7 color-coded single-channel models
- high accuracy and precision according to DIN EN ISO 8655-2
- · ergonomic designed and low weighted
- · optimized spring force for comfortable pipetting
- · robust and stable construction
- · tip ejector and piston made of stainless steel
- suitable for all established brands of pipette tips
- completely autoclavable at 121°C*
- · individual certificate and serial number
- DE-M marked according to the German calibration law

WITOPET professional

Volume	Color	Division	Volume tested	Systematic errors	Random errors	Order number
μl		μl	μl	±μl	≤ µl	
0.2 - 2.0	red	0.1	0.2 0.5 1 2	0.08 0.08 0.08 0.08	0.04 0.04 0.04 0.04	5 402 002
0.5 - 10	red	0.1	1 5 10	0.12 0.12 0.12	0.08 0.08 0.08	5 402 010
2 - 20	yellow	0.1	2 10 20	0.20 0.20 0.20	0.10 0.10 0.10	5 402 020
10 - 100	yellow	1.0	10 50 100	0.80 0.80 0.80	0.30 0.30 0.30	5 402 100
20 - 200	yellow	1.0	20 100 200	1.60 1.60 1.60	0.60 0.60 0.60	5 402 200
100 - 1000	blue	10	100 500 1000	8 8 8	3 3 3	5 402 901
1000 - 5000*	black	10	1000 2500 5000	40 40 40	15 15 15	5 402 950

Pipette support for WITOPET professional pipettes

5 406 107

^{*5000}ml filter is not autoclavable