

KOMPATIBILITÄTSTABELLE THERMOPIPETTES

OMNITIP TIPS	TIP VOLUME [µL]	CAT NUMBER	Finnpipette Digital single channel						Finnpipette Digital multichannel			Finnpipette Focus single channel					Finnpipette F1/F2/F3 single channel					Finnpipette F1/F2/F3 multichannel							
			0,2-10	0,5-10	2-20	10-100	20-200	100-1000	0,5-10	5-50	50-300	1-10	5-50	10-100	30-300	100-1000	0,2-2	1-10	2-20	10-100	20-200	30-300	100-1000	500-5000	1-10	5-50	10-100	30-300	
FILTER	0,5-10	81240	*						*			*					*								*				
	0,5-20	86240																											
	0,5-50	87240				*				*			*						*							*			
	1-100	82240				*	*						*						*	*							*		
	1-200	83240									*			*							*						*		*
	100-1000	85240						*							*							*	*						

OMNITIP TIPS	TIP VOLUME [µL]	CAT NUMBER	Finnpipette Digital single channel						Finnpipette Digital multichannel			Finnpipette Focus single channel					Finnpipette F1/F2/F3 single channel					Finnpipette F1/F2/F3 multichannel							
			0,2-10	0,5-10	2-20	10-100	20-200	100-1000	0,5-10	5-50	50-300	1-10	5-50	10-100	30-300	100-1000	0,2-2	1-10	2-20	10-100	20-200	30-300	100-1000	500-5000	1-10	5-50	10-100	30-300	
NON FILTER	0,5-10	81011	*	*					*								*								*				
		81110	*	*					*								*								*				
		81210	*	*					*								*								*				
		81220	*	*					*								*								*				
		81311	*	*					*								*								*				
		81410	*	*					*								*								*				
		81510	*	*					*								*								*				
		81710	*	*					*								*								*				
	1-200	84110			*	*	*				*			*	*					*	*	*					*	*	
		84210			*	*	*				*			*	*					*	*	*					*	*	
		84220			*	*	*				*			*	*					*	*	*					*	*	
		84311			*	*	*				*			*	*					*	*	*					*	*	
		84410			*	*	*				*			*	*					*	*	*					*	*	
		84510			*	*	*				*			*	*					*	*	*					*	*	
		84710			*	*	*				*			*	*					*	*	*					*	*	
	1-300	84110			*	*	*				*			*	*					*	*	*					*	*	
		84210			*	*	*				*			*	*					*	*	*					*	*	
		84220			*	*	*				*			*	*					*	*	*					*	*	
		84311			*	*	*				*			*	*					*	*	*					*	*	
		84410			*	*	*				*			*	*					*	*	*					*	*	
		84510			*	*	*				*			*	*					*	*	*					*	*	
		84710			*	*	*				*			*	*					*	*	*					*	*	
		85011			*	*	*				*			*	*					*	*	*					*	*	
	100-1200	85110						*							*								*						
		85210						*							*								*						
		85220						*							*								*						
		85311						*							*								*						
		85410						*							*								*						
85510							*							*								*							
85710							*							*								*							