Reax 2 inlc. universal adaptor

P/N: 541-21009-00



- Quickly clamps in 2 vessels of any size and accepts load capacities up to 1 kg
- Speed continuously adjustable from 20 to 100 rpm
- Universal adaptor allows for the use of various vessels of 50 160-mm height
- In addition, an optional adaptor is available for the use of 20 test tubes at a time

Reax 2 inlc. universal adaptor - Technical Data

Overheat protection	self-resetting
Permissible ambient conditions	5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity
Weight	5.2 kg
Protection class IEC 60529	IP 21
Timer	-
Rotation speed range	20 - 100 rpm
Power input	27 W
Dimensions (w/d/h)	510 x 180 x 235 mm
Max. Load	1 kg
Operating Mode	-
Rotation speed setting	analog
Orbit / Stroke	-
Platform size	W 230 x D 198 mm
Accessories included	universal adapter
Angle	-
Motion	overhead