





VCV 1 Vacuum control valve manual

/// Data Sheet

Allows for adjusting the set vacuum value on the manual manometer.

VC 10.400 stand not included in the scope of delivery, available as accessory (ID no. 0035001735).







designed for scientists

Technical Data

Connection diameter suction side [mm]	8	
Connection diameter pressure side [mm]	8	
Input pressure [mbar]	1 - 1050	
Pressure unit / scale	mbar	
Control range [mbar]	1 - 1050	
Resolution pressure [mbar]	50	
Measurement uncertainty [mbar]	25	
Medium temperature (gas) [°C]	5 - 40	
Venting valve	yes	
Material in contact with medium	Al2O3, PP, stainless steel	
Housing material	PP	
Fastening	stand / clamp	
Dimensions (W x H x D) [mm]	65 x 255 x 140	
Weight [kg]	0.6	
Permissible ambient temperature [°C]	0 - 40	
Permissible relative humidity [%]	80	
Protection class according to DIN EN 60529	IP 65	



