



designed for scientists



S 50 KG - HH - G 45 G - ST Dispersing tool

/// Data Sheet

Suitable for dispersing instrument T 50 digital

Ultimate fineness, emulsions 10 - 30 μm^*

Working range: 0,5 - 20 l

Stator diameter: 45 mm

Immersion depth: 70 - 260 mm

Materials in contact with medium: FFPM, SIC, AISI 316L

Suitable for solvents: yes

Max. temperature: 180 °C

Sterilization methods: all methods



designed for scientists

Technical Data

Weight [kg]	0.576
-------------	-------