

# Universal tray KM



## Description

Coated tray with drillings for individual loading with spring clamps or test tube racks.

### Loading Capacity Universal Tray KM Spring clamps (stainless steel)

| Size          | 10 ml | 25 ml | 50 ml | 100 ml | 250 ml | 500 ml | 1000 ml | 2000 ml |
|---------------|-------|-------|-------|--------|--------|--------|---------|---------|
| Pieces (max.) | 60    | 30    | 30    | 16     | 12     | 9      | 5       | 2       |

### Loading Capacity Universal Tray KM Test tube racks (with hinged foot)

Pieces max: 3

Please note:

For fastening test tube racks on a universal tray **test tube rack** and **hinged foot** are necessary.

## Details

**Shaking platform:** 356 x 300 mm

**Order No.:** 0052 030

**Accessory:**

- Spring clamps 10 ml
- Spring clamps 25 ml
- Spring clamps 50 ml
- Spring clamps 100 ml
- Spring clamps 250 ml
- Spring clamps 500 ml
- Spring clamps 1000 ml
- Spring clamps 2000 ml
- Test tube rack with 44 drillings à 14 mm Ø
- Test tube rack with 44 drillings à 16 mm Ø
- Test tube rack with 44 drillings à 18 mm Ø
- Test tube rack with 14 drillings à 30 mm Ø
- Hinged foot for test tube racks