

Universal tray KM Mini



Description

Coated tray with drillings for individual loading with spring clamps or test tube racks.

Loading Capacity Universal Tray KM Mini Spring clamps (stainless steel)

Size	10 ml	25 ml	50 ml	100 ml	250 ml	500 ml	1000 ml	2000 ml
Pieces (max.)	23	15	15	7	4	2	1	1

Loading Capacity Universal Tray KM Mini Test tube racks (with hinged foot)

Pieces max: 1

Please note:

For fastening test tube racks on a universal tray **test tube rack** and **hinged foot** are necessary.

Details

Shaking platform: 200 x 295 mm

Order No.: 0052 119

Accessory:

- Spring clamps 10 ml
- Spring clamps 25 ml
- Spring clamps 50 ml
- Spring clamps 100 ml
- Spring clamps 250 ml
- Spring clamps 500 ml
- Spring clamps 1000 ml
- Spring clamps 2000 ml
- Test tube rack with 44 drillings à 14 mm Ø
- Test tube rack with 44 drillings à 16 mm Ø
- Test tube rack with 44 drillings à 18 mm Ø
- Test tube rack with 14 drillings à 30 mm Ø
- Hinged foot for test tube racks