



## Labworldsoft® 6 Lite

### /// Data Sheet

Labworldsoft® 6 Lite laboratory software enables the monitoring, control and data recording of your IKA device. The Lite version is the entry-level model from the labworldsoft® family and ideally suited to monitor and control a single IKA device. With the help of a PC laboratory processes can be automated, measured value recordings can be defined or real-time calculations can be made.

labworldsoft® 6 Lite at a glance

- > Ideal for recording measured values & controlling one IKA device
- > 1 IKA device from the device library can be selected



designed for scientists

- > Supports requirements of FDA 21 CFR Part 11
- > Supports advanced user management (user, super user, administrator)
- > Supports labworldsoft® server
- > 1 year of support with new versions of labworldsoft® 6 Starter

#### Your benefits

- > Time savings and increased safety through automation of processes
- > Great flexibility through individual configuration of laboratory processes

#### Hardware and software requirements:

- > Operating System: Windows 10 May 2019 Update (version 1903) 64-bit or later / Windows 11 64-bit
- > .NET Framework: .NET Framework 4.8 or later
- > at least 8 GB RAM
- > at least 10 GB free hard disk space

#### Try it free of charge

You can test labworldsoft® 6 for 30 days free of charge: During this time, you have access to the entire device library with over 700 laboratory devices as well as to all functions that are required for your new digital laboratory. There is no obligation for you. After 30 days, you either choose one of our licenses or your access will expire.



designed for scientists