

XS pH Tester 1 - KIT

Code: 50014013



TECHNICAL DATA SHEET

Supplied as standard

Standards applied

Regulations	CE
-------------	----

Performance

pH measuring range	0,0 ... 14,0
pH resolution	0.1

Features

Display	LCD
IP rating	67
PH calibration points	1 ... 2 (choice of 4.01 - 7.00 - 10.01)
Recognised pH buffers	3 US pads

Feeds

Power supply	2 x AAA batteries 1.5 V
--------------	-------------------------

Dimensions and weight

Dimensions (WxDxH)	40 x 200 mm
Weight	0.105 Kg