

XS Tester pH 5 - ECOPACK

Code: 50014053



TECHNICAL DATA SHEET

Supplied as standard

Standards applied	
Regulations	CE
Performance	
Temperature accuracy	±0,2 °C
pH measuring range	-2,00 ... 16,00
temperature measurement range	0,0 ... +60,0 °C
pH resolution	0.01
Temperature Resolution	0.1 °C
Features	
Display	3-colour backlit LCD
IP rating	67
PH calibration points	1 ... 3
Recognised pH buffers	5 US pads
Feeds	
Power supply	2 x AAA batteries 1.5 V
Dimensions and weight	
Dimensions (WxDxH)	40 x 200 mm
Weight	0.11 Kg