

XS Tester pH 5 Food - KIT

Code: 50014113



TECHNICAL DATA SHEET

Supplied as standard

Standards applied

Regulations	CE
-------------	----

Performance

Temperature accuracy	±0,2 °C
pH measuring range	-2,00 ... 16,00
temperature measurement range	0,0 ... +60,0 °C
pH resolution	0.01
Temperature Resolution	0.1 °C

Features

Display	3-colour backlit LCD
IP rating	67
PH calibration points	1 ... 3
Recognised pH buffers	5 US pads

Feeds

Power supply	2 x AAA batteries 1.5 V
--------------	-------------------------

Dimensions and weight

Weight	0.14 Kg
--------	---------