

Via della Meccanica, 25, 41012 Carpi (MO) - Italy
+39 059 653274 - info@giorgiobormac.com

XS Sensor Standard DIN

Code: 32200273



TECHNICAL DATA SHEET

Performance

Measuring range pH	0 ... 14 pH
temperature measurement range	0 ... +100 °C
mV measurement range	

Features

Certificate	Yes
Series	DIN
Connection	Fixed cable / DIN
Cable	Si
cable lenght	1 m
Diaphragm	Ceramic
Electrolyte	KCl 3M
Body material	Glass
Serial number	Si
Reference system	cartridge/encapsulated

Dimensions

Body diameter	12 mm
Body lenght	120 mm