

XS Sensor Standard T DIN

Code: 32200463



TECHNICAL DATA SHEET

Performance

Measuring range pH	0 ... 14 pH
temperature measurement range	-10 ... +100 °C
mV measurement range	

Features

Series	DIN + BANANA
Connection	Fixed cable / DIN
Cable	Si
cable lenght	1 m
Diaphragm	Ceramic
Electrolyte	KCl 3M
Body material	Glass
Temperature sensor	Si / spina banana 4 mm
Reference system	cartridge/encapsulated

Dimensions

Body diameter	12 mm
Body lenght	120 mm