

XS Sensor Flow T DIN PT 1000

Code: 32200733



TECHNICAL DATA SHEET

Performance

Measuring range pH	0 ... 14 pH
temperature measurement range	-10 ... +80 °C
mV measurement range	

Features

Certificate	Yes
Series	DIN
Connection	Fixed cable / DIN
Cable	Si
cable lenght	1 m
Diaphragm	Collar
Electrolyte	KCl 3M
Body material	Glass
Serial number	Si
Temperature sensor	Si / spina banana 4 mm
Reference system	cartridge/encapsulated

Dimensions

Body diameter	12 mm
Body lenght	120 mm