

XS Sensor 2-Pore T NTC with 1m BNC Cable

Code: 32200323



TECHNICAL DATA SHEET

Performance

Measuring range pH	0 ... 14 pH
temperature measurement range	0 ... +60 °C
mV measurement range	

Features

Series	BNC + RCA
Connection	BNC + RCA
Cable	Si
cable lenght	1 m
Diaphragm	Open hole
Electrolyte	Polymer
Body material	Glass
Temperature sensor	Si / spina RCA Cinch
Reference system	Encapsulated ag

Dimensions

Body diameter	12 mm / 6 mm tip
Body lenght	95 mm / 35 mm