

XS Sensor Food T BNC

Code: 32200753



TECHNICAL DATA SHEET

Performance

Measuring range pH	0 ... 14 pH
temperature measurement range	-10 ... +100 °C
mV measurement range	

Features

Series	BNC + RCA
Connection	BNC + RCA
Cable	Si
cable lenght	1 m
Diaphragm	Ceramic
Electrolyte	Glycerin KCl
Body material	Glass
Temperature sensor	Si / spina RCA Cinch
Reference system	cartridge/encapsulated

Dimensions

Body diameter	12 mm
Body lenght	120 mm