

XS Sensor Zeta Trode S7

Code: 32200403



TECHNICAL DATA SHEET

Performance

Measuring range pH	0 ... 14 pH
temperature measurement range	0 ... +60 °C
mV measurement range	

Features

Series	S7
Connection	Screw head S7
Diaphragm	Ceramic
Electrolyte	Glycerin KCl
Body material	Glass
Reference system	Ag /AgCl

Dimensions

Body diameter	12 mm / 4 mm tip
Body lenght	190 mm / 120 mm