

XS Sensor 2-Pore F TEMP DIN (Pt1000)

Code: 32200713



TECHNICAL DATA SHEET

Performance

Measuring range pH	0 ... 14 pH
temperature measurement range	0 ... +60 °C
mV measurement range	

Features

Series	DIN + BANANA
Connection	DIN + Banana
Cable	Si
cable lenght	1 m
Diaphragm	Open hole
Electrolyte	Polymer
Body material	Plastic (ABS/POM) / Glass
Temperature sensor	Si / spina banana 4 mm
Reference system	Encapsulated ag

Dimensions

Body diameter	12 mm / 6 mm tip
Body lenght	92 mm / 50 mm