

XS Sensor 2-Pore K Bladeless

Code: 32200343



TECHNICAL DATA SHEET

Performance

Measuring range pH	0 ... 14 pH
temperature measurement range	0 ... +60 °C
mV measurement range	

Features

Series	WITHOUT BLADE
Connection	Screw head S7
Diaphragm	Open hole
Electrolyte	Polymer
Body material	Atainless Steel / Glass
Reference system	Encapsulated ag

Dimensions

Body diameter	12 mm / 6 mm tip
Body lenght	110 mm / 35 mm