

XS Sensor 2-Pore steel S7

Code: 32200653



TECHNICAL DATA SHEET

Performance

Measuring range pH	o ... 14 pH
temperature measurement range	o ... 80 °C
mV measurement range	

Features

Series	S7
Connection	Screw head S7
Diaphragm	Open hole
Electrolyte	Polymer
Body material	Atainless Steel / Glass
Reference system	Encapsulated ag

Dimensions

Body diameter	12 mm / 6 mm tip
Body lenght	105 mm / 50 mm