

XS Sensor 2-Pore steel T

NTC+1MT BNC

Code: 32200313



TECHNICAL DATA SHEET

Performance

Measuring range pH	o ... 14 pH
temperature measurement range	o ... 80 °C
mV measurement range	

Features

Series	BNC + RCA
Connection	BNC + RCA
Cable	Si
cable lenght	1 m
Diaphragm	Open hole
Electrolyte	Polymer
Body material	Atainless Steel / Glass
Temperature sensor	Si / spina RCA Cinch
Reference system	Encapsulated ag

Dimensions

Body diameter	12 mm / 6 mm tip
Body lenght	105 mm / 50 mm