

XS Sensor 2 Pore F Temp-BNC

DHS

Code: 32200143



TECHNICAL DATA SHEET

Performance

Measuring range pH	0 ... 14 pH
temperature measurement range	0 ... +60 °C
mV measurement range	

Features

Series	DIGITAL pH ELECTRODE
Connection	Fixed cable / BNC
Cable	Si
cable lenght	1 m
Diaphragm	Open hole
Electrolyte	Polymer
Body material	Plastic (ABS/POM) / Glass
Temperature sensor	Si / spina RCA Cinch
Reference system	Covered with Ag

Dimensions

Body diameter	12 mm / 6 mm tip
Body lenght	90 mm / 50 mm