

XS Sensor Electrode 201T DHS

Code: 32200103



TECHNICAL DATA SHEET

Performance

Measuring range pH	0 ... 14 pH
temperature measurement range	0 ... +60 °C
mV measurement range	

Features

Series	DIGITAL pH ELECTRODE
Connection	Fixed cable / BNC
Cable	Si
cable lenght	1 m
Diaphragm	Ceramic
Electrolyte	Gel
Body material	Plastic
Temperature sensor	Si / spina RCA Cinch
Reference system	cartridge/encapsulated

Dimensions

Body diameter	12 mm
Body lenght	120 mm