

XS Sensor Flow Temp-BNC

DHS

Code: 32200133



TECHNICAL DATA SHEET

Performance

Measuring range pH	0 ... 14 pH
temperature measurement range	-10 ... +80 °C
mV measurement range	

Features

Certificate	Yes
Series	DIGITAL pH ELECTRODE
Connection	Fixed cable / BNC
Cable	Si
cable lenght	1 m
Diaphragm	Collar
Electrolyte	KCl 3M
Body material	Glass
Serial number	Si
Temperature sensor	Si / spina RCA Cinch
Reference system	Covered with Ag

Dimensions

Body diameter	12 mm
Body lenght	120 mm