

XS Sensor Standard T DHS BNC

Code: 32200123



TECHNICAL DATA SHEET

Performance

Measuring range pH	0 ... 14 pH
--------------------	-------------

Features

Certificate	Yes
Series	DIGITAL pH ELECTRODE
Connection	Fixed cable / BNC
Cable	Si
cable lenght	1 m
Diaphragm	Ceramic
Electrolyte	KCl 3M
Body material	Glass
Serial number	Si
Temperature sensor	Si / spina RCA Cinch
Reference system	cartridge/encapsulated

Dimensions

Body diameter	12 mm
Body lenght	120 mm