

XS Sensor Standard Ag S7

Code: 32200483



TECHNICAL DATA SHEET

Performance

temperature measurement range	-10 ... +80°C
mV measurement range	-2000 ... +2000 mV

Features

Series	S7
Connection	Screw head S7
Diaphragm	Ceramic
Electrolyte	KNO ₃ 1M
Body material	Glass
Temperature sensor	-10 ... +80 °C
Reference system	Cartridge/encapsulated

Dimensions

Body diameter	12 mm
Body lenght	120 mm