

MS-D Digital magnetic ARGOlab stirrer

Code: 22002003



TECHNICAL DATA SHEET

Performance

Working humidity	80 % RH
Stirring speed	200 - 1500
Max stirring quantity	3 L

Features

IP rating	42
Speed indication	Digital via LED display
Maximum length of magnetic anchor	25
Working plate - material	Nylon+GF

Feeds

Power supply	220
Frequency	50 / 60

Dimensions and weight

Dimensions (WxDxH)	260 x 150 x 80 mm
Weight	1.8
Working plate - dimensions	135 mm