

MULTISPIN 12 - Mini

Microcentrifuge 6x2ml,

12,000rpm

Code: 40005012



TECHNICAL DATA SHEET

Standards applied

Warranty	2 years
----------	---------

Performance

Capacity max. (fixed angle rotor)	6 x 2 ml
RCF max (fixed angle rotor)	7900 xg
Noise	45 dB
Maximum speed with fixed angle rotor	12000 rpm

Features

Control system	Microprocessor and LED display
Timer	1 - 99 min hh:mm
Type of motor	Brushless induction

Feeds

Power supply	220V, 50-60 Hz
Consumption	30 W

Dimensions and weight

Dimensions (WxDxH)	176 x 196 x 122 mm
Weight	1.5 kg