

Code: 40000062



TECHNICAL DATA SHEET

Standards applied

Warranty	2 years
----------	---------

Performance

Capacity max. (fixed angle rotor)	12 x 15 ml
Deceleration ramps	BREAK function
RCF max (fixed angle rotor)	4280 xg
Noise	53 dB
Maximum speed with fixed angle rotor	5250 rpm

Features

Imbalance control	Yes
RCF display	No
Spin function	No
RCF setting	No
Programs	No
Control system	Microprocessor and LED display
Timer	Da 1 a 99 min o in continuo hh:mm
Type of motor	Brushless induction

Feeds

Power supply	220 V ±10% - 50 Hz
Consumption	250 W
Max. amperes	1 A

Dimensions and weight

Entrance height	250 mm
Dimensions (WxDxH)	320 x 440 x 290 mm
Dimensions height when open	530 mm
Weight	18,5 kg

Via della Meccanica, 25, 41012 Carpi (MO) - Italy
+39 059 653274 - info@giorgiobormac.com

Weight (with packaging, no rotor)	31 kg
-----------------------------------	-------
