

DBR 95 DIGITAL REFRACTOMETER 0-95 Brix

Code: 44000063



TECHNICAL DATA SHEET

Performance

Brix accuracy	±0,3 %
Refractive Index Accuracy	±0,0002 nD
Temperature accuracy	±0,3 °C
Temperature compensation	Automatic
Brix measuring range	0,0 ... 95,0 % Brix
Refraction index measuring range	1,3330 ... 1,5318 nD
temperature measurement range	0,0 ... 40,0 °C
Brix resolution	0.1 %
Temperature Resolution	0.1°C
Working temperature	0,0 ... 40,0 °C

Features

Feeds

Power supply	1 x lithium-ion battery 3.7 V
--------------	-------------------------------

Dimensions and weight

Dimensions (WxDxH)	180 x 100 x 55
Weight	0.41 kg