

DBR 55 DIGITAL REFRACTOMETER 0-55 Brix

Code: 44000073



TECHNICAL DATA SHEET

Performance

Brix accuracy	± 0,2 %
Temperature accuracy	±1,0 °C
Temperature compensation	Yes, 0 ... + 60 °C automatic
Brix measuring range	0,0 ... 55,0 % Brix
temperature measurement range	0,0 ... 70,0 °C
Brix resolution	0.1 %
Temperature Resolution	0,1 °C
Working temperature	10,0 ... 60,00 °C

Features

Feeds

Dimensions and weight

Dimensions (WxDxH)	58 x 36 x 122 mm
Weight	0.15 kg