

MOD.110 ABBE REFRACTION Meter

Code: 43000133



TECHNICAL DATA SHEET

Performance

Brix accuracy	$\pm 0,25 \%$
Refractive Index Accuracy	$\pm 0.0002 \text{ nD}$
Brix measuring range	0,00 ... 95,00 % Brix
Refraction index measuring range	1.300 ... 1.700 nD
temperature measurement range	0 ... 70 °C

Features

Feeds

Dimensions and weight

Dimensions (WxDxH)	100 x 200 x 240 mm
Weight	2.6 kg