

M2-D PRO ARGOLab Digital heating stirrer

Code: 22002063



TECHNICAL DATA SHEET

Performance

Heating power	500
Over Temperature Safety Device	50
Working temperature	5 ... 40
Heating Temperature Range	Ta ... 340 °C
Working humidity	80 % RH
Stirring speed	100 ... 1500
Max stirring quantity	20 L

Features

IP rating	42
Temperature indication	Digital
Speed indication	LCD Display
Maximum length of magnetic anchor	80
Working plate - material	Ceramic-coated stainless steel

Feeds

Frequency	50 / 60
-----------	---------

Dimensions and weight

Dimensions (WxDxH)	280 x 160 x 65 mm
Weight	2.8
Working plate - dimensions	Ø135 mm